

DIGITAL BABY LENGTH MEASURING DEVICE



DETECTO
Model DLM

**Readouts in feet/inches,
inches, and centimeters.**

- May be used stand-alone on an exam table.
- (2) AAA batteries included.

***Easy-to-read, bold LCD with convenient units switching.
HOLD feature retains measurement for active babies.***

DETECTO DIGITAL BABY LENGTH MEASURING



Detecto's model DLM Digital Baby Length Measuring Device provides lightning-fast readouts. The lightweight design and battery power make it ideal for easy transport from one patient room to another.

DLM

DIGITAL SPEED + ACCURACY

Detecto model DLM



CONVENIENT UNITS SWITCHING

Choose from ft/in, in, or cm.

BATTERY POWER

Detecto's DLM comes with two 1.5V AAA batteries making it completely portable for easy mobility.



DLM SPECIFICATIONS

MEASUREMENT RANGE:	13¾ in - 31½ in / 35 cm - 80 cm
GRADUATIONS:	1/16 in / 1 mm
OVERALL DIMENSIONS:	11.4 in W x 2.7 in H x 24.4 in L / 29 cm W x 7 cm H x 62 cm L
POWER:	(2) 1.5V AAA batteries (included)
FUNCTION KEYS:	RESET, HOLD, UNIT, ON/OFF
SHIPPING WEIGHT:	2 lb / 907 g
TEMPERATURE RANGE:	+50°F up to 104°F / +10°C up to +40°C

- ▶ **Rubber Feet Ensure Stable Footing While Measuring**
- ▶ **Smooth, Contoured Paddles are Gentle to Infant Touch**

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

Sold By:

**DETECTO**
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com