

# Suspended Baby Scale



- Ultraportable
- Accurate
- No Power Required

*Detecto's model  
HS25KGP  
Suspended  
Baby Scale*

***Suspended Sling  
Seat Included***

Ideal for doctors and nurses working in remote locations, Detecto's model HS25KGP Baby Scale is lightweight and portable without requiring any battery or AC power. The waterproof vinyl weighing trousers are attached to a heavy-duty, zinc-plated S-hook providing a comfortable suspended sling seat for weighing babies. The 20 cm diameter, shatterproof dial features a white face with black graduations and numerals for high-contrast readouts.

Model	Capacity	Shipping Weight
HS25KGP	25 kg x 100 g	5 lb / 2 kg