

## MB130 | DIGITAL PEDIATRIC SCALE WITH WEIGHING TRAY



DETECTO model MB130 scale with optional DLM digital length measuring

**COMFORTABLE  
CONTOURED  
WEIGHING TRAY  
CRADLES BABIES SECURELY**



### DLM

DETECTO's optional DLM digital length measuring may be added to the MB130 scale providing a completely digital measuring station for pediatrics.



### MB-CASE

Optional water-resistant carrying case with adjustable shoulder strap available for the MB130 scale made of robust nylon (model MB-CASE).

**Babies feel at ease and safe in the extra-large weighing tray**

- 44 lb / 20 kg Capacity and 0.01 lb / 5 g Accuracy
- Portable and Lightweight for Mobile Medical Use
- Easy to Operate with Large Readouts
- Integral Analog Measuring Tape
- Optional Digital Length Measuring Available



## FROM NEONATAL TO PEDIATRIC WEIGHING ACCURACY

DETECTO's model MB130 digital baby scale features an easy-to-read one-inch-high LCD with pounds and kilograms readouts for neonatal to pediatric measurement. Power is provided by either 9-volt battery (not included) or AC adapter (included). The auto shut-off after two minutes helps conserve battery life.

The convenient HOLD buttons on the scale (and also on the optional digital length measuring rod) retain the infant's measurements on the display even after they are taken off the scale, so the patient can be attended to first before recording measurements.

The easy-to-clean plastic tray is ideal for sterile environments and features a roomy space that is 22-inch-wide by 11-inch-deep (56-cm-wide by 29-cm-deep) for babies to lay comfortably. The surface material is kind to the skin and is easy to disinfect.

The convenient Zero function allows the user to automatically subtract the weight of a blanket, toy, or diaper from the infant's weight measurement. Weighing only 5 lbs, the MB130 is handy and lightweight for transporting to different locations. The optional carrying case eases portability and provides safe storage when not in use.



## REMOVABLE TRAY FOR STANDING TODDLER WEIGHING

Multi-functionality is a key feature of DETECTO's MB130 scale. The baby tray can be quickly removed which automatically converts the scale into a low-profile step-on toddler scale with a capacity of up to 44 pounds / 20 kilograms.

# DETECTO's MB130 PEDIATRIC SCALE

Easy to Operate with Large Tray  
and Easy-to-read Display

## DIGITAL LENGTH MEASUREMENT ACCURACY

Weigh and measure in one step



DETECTO's optional DLM digital length measuring rod can be easily mounted to your MB130 scale providing a completely digital weight and height measurement station for pediatrics. The DLM provides readouts in feet/inches, inches, and centimeters on a bold, one-inch-high LCD. The HOLD button retains measurements for active babies.

The DLM offers you digital speed and accuracy and it comes with two 1.5V AAA batteries making it completely portable. The smooth, contoured paddles are gentle to infant touch and measurements range from 13-3/4 inch – 31-1/2 inch / 35 cm – 80 cm with 1/16 inch / 1 mm graduations.



Optional Carrying Case, Digital Length  
Measuring and Baby Scale Cart Also Available

## CARRYING CASE (OPTIONAL)



DETECTO offers an optional carrying case designed to comfortably transport the MB130 baby scale. The MB-CASE features a high-quality robust nylon construction that is water and dust resistant making it ideal for mobile healthcare environments.

## MOBILE BABY SCALE CART



DETECTO's model RC33669WHT baby scale cart perfectly accommodates the MB130 scale and digital length measuring rod for a mobile pediatric measurement station throughout your clinical facility. The steel rolling cart with white powder paint features a push handle and two parking locks plus one steering lock for optimal patient safety. The top shelf is made of ABS molded plastic and will securely grip the MB130 scale's four rubber feet without sliding for patient stability. Five drawers provide ample storage for pediatric care supplies (two 3-in-high, two 6-in-high, and one 9-in-high).

# ECONOMIC, PORTABLE PEDIATRIC MEASUREMENT



When fitted with the optional DLM measuring rod, the MB130 baby scale becomes an all-digital, time-saving measuring station for weight and height.



<b>Capacity</b>	0 - 22 lb x 0.01 lb / 22 - 44 lb x 0.02 lb 0 - 10 kg x 5 g / 10 - 20 kg x 10 g	
<b>Units of Measure</b>	lb / kg	
<b>Function Keys</b>	ON/ZERO/OFF, HOLD, UNIT	
<b>Stabilization Time</b>	1-2 seconds	
<b>Operating Temperature</b>	40° F - 105° F / 5° C - 40° C	
<b>Power</b>	(1) 9V battery (not included) or AC adapter (included)	
<b>Display</b>	1 in / 25 mm high, 5 digit LCD	
<b>Measuring Tape</b>	Up to 21-9/16 in / 55 cm with 1/16 in / 1 mm increments	
<b>Dimensions</b>	<b>Scale Base Only</b>	12.2 in W x 12.6 in D x 1.3 in H / 31 cm W x 32 cm D x 3.3 cm H
	<b>Scale with Tray</b>	22 in W x 13.1 in D x 4 in H / 55.9 cm W x 33.3 cm D x 10.2 cm H
	<b>Scale with Tray and DLM</b>	24.6 in W (expands to 34.4 in W) x 16 in D x 6.3 in H / 62.5 cm W (expands to 87.4 cm W) x 40.6 cm D x 16 cm H
	<b>Weight Tray</b>	22 in W x 11.3 in D x 3.7 in H / 55.9 cm W x 28.7 cm D x 9.4 cm H
<b>Net Weight</b>	5.4 lb / 2.4 kg	
<b>Shipping Weight</b>	8 lb / 3.6 kg	
<b>UPC Code</b>	809161184402	



When purchasing a MB130 scale, consider adding the optional MB-CASE for site-to-site transport or the Rescue brand baby scale cart for rolling mobility within a clinical facility.

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.



203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

**Mobile:** [m.detecto.com](http://m.detecto.com)



Scan this QR code on your smartphone for more product info.

**SOLD BY:**