



Eye-Level Mechanical Physician Scales

Mild Steel with Electrostatic Powder Paint Finish for Optimum Long-Lasting Quality

**Now available
in 500 lb / 220 kg
capacities!**

*Detecto
model 448
featuring a
height rod,
handpost,
and wheels*



**HEAVY-DUTY
SOLID, STABLE BASE**
*for patient assurance
while weighing.*

WEIGH UP TO 500 LB / 220 KG CAPACITY

Detecto's innovatively-designed scales feature dual-readings on either side of the weigh beam making weighing simple and fast. Pounds, kilograms, or pounds/kilograms are all available on the aluminum easy-to-read beam. The dual-reading height rod measures in inches or centimeters at a glance, and the platform cover can be easily removed for cleaning.



**Specify
Detecto**

**It's the
name that
means
scales in the
health care
industry.**

The dual-reading, die-cast aluminum on black insert weighbeam is easy to read and makes weighing simple and fast.



HANDPOST

An important aid for the elderly or weak individual. Models with this feature listed below.



HEIGHT ROD

Measured in inches and centimeters: 30 - 78 inches and 76 - 200 centimeters



WHEELS

Two rear wheels allow easy transport. The light weight of the Detecto physician scales combined with optional wheels make transport simple.

Eye-Level Mechanical Physician Scales



MODEL	CAPACITY		HEIGHT ROD 30 - 78 in 76 - 200 cm	WHEELS	HANDPOST	SHIPPING WEIGHT
	Avoirdupois	Metric				
439	400 lb x 4 oz		•			38 lb/17.2 kg
2391		180 kg x 100 g	•			37 lb/16.7 kg
437	400 lb x 4 oz					35 lb/15.8 kg
2371		180 kg x 100 g				35 lb/15.8 kg
438	400 lb x 4 oz		•	•		39 lb/17.6 kg
2381		180 kg x 100 g	•	•		39 lb/17.6 kg
448	400 lb x 4 oz		•	•	•	43 lb/19.5 kg
449	400 lb x 4 oz		•		•	43 lb/19.5 kg
2491		180 kg x 100 g	•		•	43 lb/19.5 kg
339	400 lb x 4 oz	175 kg x 100 g	•			38 lb/17.2 kg
337	400 lb x 4 oz	175 kg x 100 g				35 lb/15.8 kg
338	400 lb x 4 oz	175 kg x 100 g	•	•		39 lb/17.6 kg
349	400 lb x 4 oz	175 kg x 100 g	•		•	43 lb/19.5 kg

Platform Size: 10.5"W x 14.5"D/27 cm W x 37 cm D

Scale Height: 59"/150 cm

Optional Counterweights for Increased Capacities

- CW-9** Counterweight for use with new 437, 438, 439, 448 and 449 scales (Increases weight range from 400 lb to 500 lb)
- CW-80X** Counterweight for use with new 2371, 2381, 2391 and 2491 scales (Increases weight range from 180 kg to 220 kg)



Cardinal Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



Sold By:



203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com