



# Mariner® WPS12DT

## Wet Diaper Scale



**MEETING  
NEONATAL  
FACILITY  
REQUIREMENTS,  
RIGHT  
DOWN TO THE  
GRAM**

- 192 oz x 0.05 oz / 5,500 g x 1 g
- Tare Button to Remove Dry Diaper Weight and Measure Fluid Output
- 80-hour Rechargeable Battery Pack with Automatic Shut Off Feature

- 5/8-in Deep Stainless Steel Removable Tray Aids Sample Retention and Allows Easy Clean Up
- 1.2-in/30-mm High Backlit LCD Easily Readable in Dimly-Lit Nurseries
- Scale is Completely Submersible Underwater

# DEPENDABLE ACCURACY

for  
**CLINICAL  
SETTINGS**



DETECTO's Mariner® WPS12DT wet diaper scale meets all the hygiene demands of any active medical environment.

## EASY-TO-USE, PORTABLE, AND CONVENIENT

The backlit blue LCD display offers large 1-in/25-mm high bold digits, readily visible in even the dimmest nursery conditions. The raised buttons allow for tactile functionality and easy operation. The UNIT key toggles through six different weight capacity display options in both decimals and fractions. The tare function conveniently allows the weight of a dry diaper to be removed from the total weight of a soiled diaper. The convenient 80-hour rechargeable battery pack with automatic shut-off feature conserves energy and allows the scale to be moved where it is needed the most.



Removable stainless steel tray makes sanitation effortless.

## HYGIENIC AND HEAVY-DUTY

The watertight AC plug and stainless steel construction allow for easy clean up, prevent corrosion, and reduce the risk of bacteria harboring in nooks and crannies. The 5/8-in deep removable stainless steel tray with raised edges holds items neatly on the weighing platform and makes clean up effortless. Both of these features are a necessity for messy nursery applications.

Backlit LCD is easily visible, even in darkened rooms so babies won't be disturbed.



Model	WPS12DT Wet Diaper Scales		Dimensions	Scale	7.5 in W x 9.8 in D x 2.3 in H / 190.5 mm W x 247.5 mm D x 57.0 mm H
Capacity x Increments	12 lb x 0.1 oz	192 oz x 0.05 oz	Tray	9.0 in W x 6.0 in D x 1.1 in H / 228.6 mm W x 152.4 mm D x 27.9 mm H	
	12 lb x 0.005 lb	192 oz x 1/8 oz			
	12 lb x 1/8 oz	5,500 g x 1 g			
Zero	Established on power up routine and maintained by auto-zero circuitry		Tare	100% of full scale capacity	
Display / Digits	1.2-in / 30-mm high backlit blue LCD, (5) 1-in / 25-mm high black LCD digits		Keypad Buttons	ON/OFF, UNITS/Backlight, SET, TARE	
Feet	6 rubberized feet for traction		Enclosure Rating	IP67 submersible	
Product Weight	4.9 lb / 2.2 kg		Operation Temperature	40° to 105°F (5° to 40°C)	
Shipping Weight	7.1 lb / 3.2 kg		Storage Temperature	40° to 140°F (5° to 60°C)	
Power	AC Adapter (included): 120 VAC input 9V, 100mA output - Battery (included): Ni-MH rechargeable 4.8 VDC 1800 mA		Battery Life	80 hours of continuous use with a fully recharged battery pack	
				Please Note: The WPS12DT Mariner® is not commercial dishwasher safe.	
UPC Code	809161201901		Certification	NSF Approved	

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

**DETECTO**  
*www.Detecto.com*

203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

Mobile: [m.detecto.com](http://m.detecto.com)

**SOLD BY:**